

LAB BENCH TABLES

Looking for lab tables and workstations for schools, hospitals, research centers, etc.?

We can help to design and customize all the laboratory needs. Our lab table is **easy installation**, **heavy duty**, **affordable**, and certificated with **ISO**, **CE**. Best of all, our laboratory workstations and tables are available in a variety of types to meet different applications.



ZH has been providing laboratory casework for ten years. Now, we are the biggest lab casework manufacturer in Guangzhou, China. Our products have sold to dozens of countries, and we have established long-term cooperation with many big customers; They are all satisfied with our quality and service. If you are looking for a reliable supplier, Nobody is more suitable than us.

CHECK ZH'S LAB TABLES AND CASEWORKS:



[Steel Lab Table](#)

Phone: +86-13302388256

Email: info@zhlabsfurniture.com

Address: No.3068, Liangsha Road, Chendong Village, Zhongluotan Town, Baiyun District, Guangzhou City, Guangdong Province, China



[Phenolic Resin Lab Table](#)



[Wood Lab Table](#)

FEATURES FOR ALL LAB BENCH TABLES

- Heavy duty-can load 300kgs/m²;

Only high-quality industrial material will be used, and our quality control department will follow it from raw material until packing.

Phone: +86-13302388256

Email: info@zhlabsfurniture.com

Address: No.3068, Liangsha Road, Chendong Village, Zhongluotan Town, Baiyun District, Guangzhou City, Guangdong Province, China

- Easy installation;

The screw holes are prepared, and snap-fit for the knee space board, the base cabinet is installed in our factory, all you need to do is connect it with screws, so the installation is very easier;

- Removable back panel for wires, pipes, etc., which can provide easy maintenance;
- Famous brand hinges, soft and silent; 110 degrees or 175 degrees; more than 50,000 times, strong anti-rust and corrosion resistance, reaching the international hardware industry standard.
- Drawer rails are silent, 45kg loading capacity. Soft-closing or heavy-duty design is available.
- Fast shipment-most of the modular lab casework is in stock.
- Provided with a leveling system that allows them to adjust the precise height of the furniture, or with heavy-duty wheels to move anywhere.

LAB CASEWORK OPTIONS

- Modular steel casework
- Stainless steel casework;
- Phenolic resin casework;
- Laboratory wood casework;
- We have experienced technicians and skilled workers and provide you with high-quality and efficient products. Email us directly for a fast and free quote today.